



Amyloid β A4, N-Terminus

clone 19H5

Order No.: 0064-100/bA4N-19H5

Size (μg) 100 Lot No.: 0064S



www.nanotools.de

orders & support:

email: info@nanotools.de phone: +49-7641-455 670 fax: +49-7641-455 671

02/020307F

Isotype	Species Reactivity	Applications	Mol. Weight	Ref.Cell Line	Epitope	Immunogen
IgG1, kappa	human	ELISA, WB		none	N-Terminus of Amyloid βA4; crossreacts with APP	N-terminal peptide conjugated to KLH

Background and Specificity:

The beta-amyloid peptide (beta A4), proteolytically released from the amyloid precursor protein (APP), is the principal component of senile plaques in Alzheimer's disease. Cleavage of APP by alpha-secretase or alternatively by beta-secretase leads to generation and extracellular release of soluble APP peptides, S-APP-alpha and S-APP-beta, respectively, and the retention of corresponding membrane-anchored C-terminal fragments, C83 and C99. Subsequent processing of C83 by gamma-secretase yields P3 peptides. This is the major secretory pathway and is nonamyloidogenic. Alternatively, presenilin/nicastrin-mediated gamma-secretase processing of C99 releases the amyloid beta proteins, amyloid-beta 40 (Abeta40) and amyloid-beta 42 (Abeta42), major components of amyloid plaques, and the cytotoxic C-terminal fragments, gamma-CTF(50), gamma-CTF(57) and gamma-CTF(59).

Mab βA4N-19H5 recognizes the free N-Terminus of the β A4 polypeptide. The antibody crossreacts with APP.

Purification: The antibody was purified from serum-free cell culture

supernatant by subsequent thiophilic adsorption and size

exclusion chromatography.

Formulation: Iyophilized from 1 ml 2 x PBS / 0.1 % Na-azide / PEG and

Sucrose.

Reconstitution: Reconstitute with 1 ml H₂O (15 min, RT).

Stability: For long-term storage, freeze lyophilizate upon arrival (-20^C).

Upon reconstitution, aliquote and freeze in liquid nitrogen; reconstituted antibody can be stored frozen at -80°C up to 1 year. Thaw aliquots at 37°C. Thawed aliquots may be stored at 4°C up to

3 months.

Avoid repeated freeze / thaw cycles.

Positive Control: none

Immunoblotting: 1 μg/ml for HRPO/ECL detection

Recommended blocking buffer: Casein/Tween 20 based blocking and blot incubation buffer, e.g. nanoTools product

#3031-500/CPPT or #3031-3000/CPPT

Immunoprecipitation: ND Immunocytochemistry: ND

ELISA: use at $0.05 \mu g/ml$

All products are supplied for research and investigational use only. Not for use in humans or laboratory animals.

Related Products

mab to β A4, N-Terminus #0084-100/bA4N-19H11

mab to βA4, N-Terminus #0195-100/bA4N-7F4

mab to βA4, N-Terminus #0196-100/bA4N-7F9

mab to βA4, N-Terminus

mab to βA4 (1-40), C-Terminus

#0060-100/bA4(40)-5C3

mab to βA4 (1-40/42), C-Terminus

#0062-100/bA4(40/42)-9F1

mab to βA4 (1-42), C-Terminus

#0061-100/bA4(42)-8G7

mab to βA4 (1-43), C-Terminus #0095-100/bA4(43)-6G12

1 2



Immunoblot Analysis

Amyloid beta A4 peptides (lane 1: bA4(1-40); lane 2: bA4 (1-42); lane 3: bA4 (1-43)) were applied on SDS-PAGE and transferred to a PVDF membrane. The immunoblot was probed with 2µg/ml mab bA4N-19H5 for 1h at 15-22°C and developed by ECL (exposure time: 30 sec).